

pcLink for Windows Reference Guide

Function Keys (function/shift function)

- F1**-Exit/Print
- F2**-Up arrow
- F3**-Down arrow
- F4**-Select
- F5**-Cut/Paste notes & schedule blocking
- F6**-Search
- F7**-Merge/Stikum
- F8**-Save screen and go to Main Menu/Calculator
- F9**-Schedule/Today's Patients
- F10**-Save screen and go to Intercom/Phone list
- F11**-Restore saved screen or return to Main Menu/Calculator

Function	Keyboard Key(s)
Select	Ctrl+P
Exit	Ctrl+[
Add Record	Ctrl+B
Remove Record	Ctrl +Z
Character, Delete Current	Ctrl +D
Character, Delete Previous	Backspace
Character Insert	Ctrl +F
Page, Next	Ctrl +T or Page Down
Page, Previous	Ctrl +R or Page Up
Page, First	Home
Page, Last	End
Restore/Main Menu	Ctrl +G
Delete Beyond Cursor	Ctrl +Y















Note: to achieve the control key functions, press and hold the **Ctrl** key and then press the corresponding key once.

pcLink Menu Bar

The following pcLink-specific functions are available on the pcLink Menu Bar:

File	Merge	Full Screen
	Print	Screen Shot
	Printer Setup	Disconnect
	Customize	
Edit	Stikum	Paste
	Cut	Calculator
Phone	Lookup	
	Intercom	
Schedule	Today's Patients	Search for Openings
	Today's Schedule	Other
Save	Save screen and go to the Main Menu.	
Restore	Restore saved screen or return to the Main Menu.	
Exit	Exit back one screen.	
Go...	User-defined quick links to integrated applications or documents.	
Help	Datacon Dental System	
	Connection	
	About	

Toolbar Buttons

	Exit back one screen.
	Launch the Merge function from the current record.
	Print. Resulting options are context-sensitive.
	Attach a stikum note to a screen.
	Search. Resulting options are context-sensitive.
	Look up patient phone numbers.
	Open the Intercom.
	Open the list of today's patients.
	Open today's schedule.
	Open the Appointment Search screen.
	Open the Scheduling screen.
	Save screen and go to the Main Menu.
	Restore saved screen or return to the Main Menu.
	Open the Calculator function.

Quit pcLink

To quit pcLink, log out of the Dental System to the System Menu and then click File>Disconnect on the pcLink Menu Bar.