


# Charting Overview

When you first access the charting area, the first page that is displayed is the Overview.



This page is intended to provide a quick overview with relevant information about the patient. The tabs across the top provide access to other areas of the charting system. There is an area reserved on this page that can be customized to display additional information about the patient that is important in your office. This might include account balances, family appointments, insurance, referral source, or other information.

You may click on any tooth to display details of oral exams or treatments that have been applied to the tooth. The buttons on the lower left side of this page provide access to historical views of the dentition, earlier exams, and an enlarged view of the teeth.

The History button on the Overview page allows you to display the dentition at any point in time. The initial menu options allow you to display the patient's teeth as they were at the first visit, as they are today, or at any date in the patient's ledger.

The screenshot shows the Datacon Dental Group software interface. At the top, there is a menu bar with various icons for navigation and actions. Below the menu bar, the patient's name "Joe Sample - 6/22/2014" is displayed. The main interface is divided into several sections:


- Navigation Tabs:** Overview, Medical History, Oral Exam, Perio Exam, Treatment, Planning, Referral, RX, Lab, Chart Notes.
- Dental Chart:** A detailed view of the patient's teeth, showing the upper and lower arches. The teeth are color-coded to indicate different conditions or treatments.
- History Table:** A table showing the patient's dental history. The table has columns for Date, Prv, Description, Tooth, and Surface. The first two rows are highlighted in yellow.
- Buttons:** Overview, History, Exams, Notes, Show, First Visit, Today, Development.
- Patient Information:** Date: 6/22/14, Patient Age: 34 years 1 month.
- Photo:** A small portrait of the patient, Joe Sample.

Date	Prv	Description	Tooth	Surface
6/17/04	AAA	Recall Examination		
6/17/04	AAA	Prophylaxis Adult		

If you select one or more teeth by clicking or clicking and dragging, the ledger will be updated to only display transactions for the selected teeth.

You may click on the Development button to access an alternate method of viewing the patient's history.


The Development menu under History on the Overview page provides a way to show the patient's dentition at any point in time.



The view shown above is the patient's dentition at the age of 5 displayed by clicking the 5 Years button. There are buttons for 6 Months, 2 Years, and 5 Years. You may hold down the Earlier or Later buttons to watch the changes over time. The Animate button allows you to enter a range of ages and sit back to watch the progress. We expect these options to be a great educational tool when parents are questioning what is going on with their children.

Other options on the Overview page include buttons to show the list of oral exams and perio exams, access to patient notes, and a Show button to display a large view.

The Show button, which is available on many of the charting screens, displays a larger view of the teeth that includes facial, occlusal, and lingual views of the teeth. You may draw on this and add notes at the bottom.



The edited image may then be printed on paper or to PDF to be saved in the patient folder.

**Questions?** Contact the staff at Datacon for help.  
www.datacondental.com  
Phone: (800) 773-7878  
Email: [info@datacondental.com](mailto:info@datacondental.com)