


Perio Charting

Introduction

Datacon's Periodontal Charting is a great way to present periodontal treatment plans to your patients, as well as an excellent teaching tool. Offices with workstations in the operatories are able to display the patient's mouth image. By looking at the screen, patients are able to see the potential for bone loss from not flossing. The tooth chart is graphically accurate down to the millimeter.

Getting Started:

- **Click on Today's Patients**  to get to *Today's Patients*
- **Click** the patient you are working with
- **Click** *Clinical Exam*
- **Click** down arrow for Perio Exam selection

Or

- **Click Patients & Accounts** from the **Main Menu**
- **Type** the patient's last name, first name, press ENTER
- Move the Tab at top to **Patient** and click on the patient
- **Click** Patient chart>Clinical records
- Answer "Y" to "A tooth chart is not on file for this patient. Do you want to add one?"
- **Click** the **Add** Button and answer "Y" to "Is this a perio exam?"
- If a Tooth Chart already exists and you are adding a new perio exam, **Click** the **Add** button, the system will default to today's date. **Click** the down arrow and select **Perio Examination**.

User's Charting Script – Perio Charting Setup

Each user on the system can define the order in which they enter their perio exam information. The charting script allows you to input the tooth data in the same sequence that it is done in the operatory. Without a charting script, all perio tooth fields display “tooth by tooth.” The beauty of entering exactly what you want in this area is that you can simply use your insert key if a tooth recession is present. This saves you from having to enter this information for every tooth if you were charting 1-16 lingual pockets, 1-16 lingual recession.

- **Click Users** from the *Main Menu*
- **Click** on your name
- **Click Charting Script**





FIGURE 1: USER PERIO CHARTING SCRIPT

- Type in your script and click **EXIT**


- The Charting Script is not needed to enter perio charting information. This is set up only if you want to streamline periodontal data entry. If another user charts exactly the same way, just copy that users charting script by hitting F5 to copy and Shift F5 to paste into the new script. If you make changes, your new charting script automatically appears as soon as you enter the patient's perio exam data.

Opening Perio Charting for the Patient

- **Click Today's Patients**  to get to *Today's Patients*
- **Click** on the patient you are working with
- **Click Clinical Exam**
- **Click Perio Exam**

OR

- **Click Patients & Accounts** from the **Main Menu**
- **Type** the patient's last name, first name, press ENTER
- Move the Tab at top to **Patient** and click on the patient
- **Click Patient chart>Clinical records**
- Answer "Y" to "A tooth chart is not on file for this patient. Do you want to add one?"
- **Click the Add Button** and answer "Y" to "Is this a perio exam?"



The screenshot shows the 'Perio Exam Information Screen' with the following components:

- Top Left:** Patient name 'Troy Peri' and date '1/07/13'. Below is a circular tooth chart with numbers 1-32. A blue square is in the center.
- Top Right:** 'Exam Record' section with a list of actions:
 - Exam Date: 1/07/13
 - Tooth Status as of 1/07/13 (Missing, Crown, etc.)
 - Teeth Examined: All
 - Enter Perio Exam Data
 - Readings Entered for 2,4-15,18-31
 - Review Perio Exam Data
 - Compare to Previous Readings
 - Print Exam Report
 - Display Image
 - Exam Notes..Probe check : ____ Full ____partial
 - Additional Exam Data
- Bottom:** A detailed tooth chart showing the upper and lower dental arches with yellow highlights on specific teeth.
- Bottom Left:** A '< Previous' button.

FIGURE 2: PERIO EXAM INFORMATION SCREEN

Preparing for Perio Exam

You can manually enter this information, or go directly into the charting area and chart existing restorations. If you have done work on the patient, it will automatically be updated.

- **Click *Tooth Status*.** You will be asked if the patient's tooth status has changed since the last exam report date. Answer "Yes" to this question. You have a choice of the following for each tooth:
 - Normal
 - Missing
 - Crown
 - Pontic
 - Implant
 - Implant and Crown
 - Impacted
 - Partial Fracture
 - Full Fracture

When you finish updating the *Tooth Status*, **Escape** and **Click Enter perio exam data**

- ***Teeth Examined*** default is ALL.
 - If only a portion of the mouth is to be examined, you can enter the tooth numbers or mouth quadrant. For example, LR is Lower Right, UL is Upper Left, 1-16 is teeth 1 through 16.
 - If you enter this information when you start the perio exam, the system will only allow you to enter perio data for that criteria.

Entering Perio Tooth Data

- **SELECT *Enter perio exam data***
 - The screen that appears depends if and how you set up your User Charting Script.




FIGURE 3: EXAMPLE OF NO CHARTING SCRIPT, WHICH ALLOWS YOU TO ENTER ALL CRITERIA FOR THE TEETH

- **Lingual Pockets:** Type the number for each pocket, type a period (.) or a plus (+) if there is bleeding, and press **ENTER**
 - Type 00, 01, or 2 through 15 for a total of three readings per surface. When you press **ENTER**, the system will go to the next tooth.
 - If a tooth starts to bleed after you have moved on, simply click on that tooth, and type . or + to designate bleeding.
- **Lingual Recession:** Type a number from -6 to 0 to 9, press **ENTER**
- **Buccal Pockets:** Type the number for each pocket, type a period (.) or a plus (+) if there is bleeding, and press **ENTER**
- **Buccal Recession:** Type a number from -6 to 0 to 9, press **ENTER**
- **Mobility:** Type 1 to 3, type a period (.) after the number to indicate a plus, and press **ENTER**
- **Morphology:** Click the down arrow to view options, and press **ENTER**
 - OK is the default.
 - Minimal gingival attachment.
 - Inadequate gingival attachment.
- **Furcation:** Type 1, 2 or 3 and press **ENTER**, or press **ENTER** if there is no furcation.
 - Furcation choices are Class I, II, or III. Furcations can only be entered on teeth that have multiple roots. Indicate whether the furcation is on the mesial, buccal or distal root.
 - Furcations are represented on the tooth chart by triangular symbols that print on the affected root.

- **Comments:** These comments are user defined for the practice. **Click** the down arrow to find the comment you want and press **ENTER**.
 - Please see **System Settings for Perio Charting** section of this chapter for this setup.
- **Prognosis:** **Click** the down arrow to find the one you want and press **ENTER**
 - *None, Good, Fair, Guarded, Very guarded, Poor, Very poor, Hopeless*
- **Suppuration:** **Click** the down arrow to choose the option and press **ENTER**.
 - The option is either *Present* or blank
 - Suppuration is a sign of infection
- **Exam notes:** This is a free form notes area for typing a narrative regarding this exam. Press **ENTER** when finished and the exam criteria will be ready for you to fill out on the next tooth.
- **EXIT** when you are finished with all the teeth and the *Exam Screen* will display.



Review Perio Exam Data

You can review and make changes to the Perio Exam with this option.

- **Click Review perio exam data**
- **Click** on a tooth number to make changes

Patient: Troy Peri Date: 4/16/13

Comment Prognosis Mobility Furcation Recession F/L Facial Lingual	01B1 3 / 5'3'4 6 6 6	1 1 1 1 1 1	3 / 4'3'3' 5 4 4	3 / 4'3'3' 6'3'4	3 /1 5 4 4 3 3 4	3 / 4 3 4 4 3 3	4 3 3 4 3 3 4 3 4	3 2 3 4 3 4 4 3 4	4 3 4 4 4 4 4'3'4'	2 / 4 3 6 5'4'4'	4 3 4 4'3'5' 4 4 5	01 /3 5 4 4 4 4 5				
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
	32	31	30	29	28	27	26	25	24	23	22	21	20	19	18	17
Lingual Facial Recession F/L Furcation Mobility Prognosis Comment	4 4 6 4 4 4 1 \	4 4 5 5 4 4 1 \	5 4 7 4 3 5 1 \	8'4'4' 4'3'4' 4 \	4 4 6 4'4'6' 4 \	4 4 4 4 4 5 2 \3	4 4 4 4 4 4 2 \3	4 4 5 4 3 4 2 \3	5 4 5 4 4 4 3 \2	4 4 4 4 3 4 3 \2	6 4 5 4 4 4 4 \	6 4 4 4 4 6 3 \	5 4 6 6 5 4 B111 B1			

< Previous

Exit

FIGURE 4: REVIEW PERIO EXAM SCREEN

Compare to Previous Readings

This option compares readings from the current chart to readings from prior visits. A change in the pocket readings (plus or minus) is highlighted. On a patient's first chart, there will be no previous readings for comparison, and any pocket higher than a two will be highlighted. You can print this screen by clicking the **PRINT** key. You can also review previous readings by pressing the page up key (or clicking on the "< Previous" tab. The numbers highlighted in red, refer to pockets that have increased from the last time you saw the patient, the numbers in green have decreased

Troy Peri															
4/30/10	4 5 6	4 5 6	2 2 3	5 6 6	4 5 3	9 5 4	4 4 4	5 8 6		5 5 6	4 5 4	4 3 4		4 3 4	
5/13/10	4 5 6	4 5 6	2 2 3	5 6 6	4 5 3	9 5 4	4 4 4	5 8 6		5 5 6	4 5 4	4 3 4		4 3 4	
5/14/10	3 4 3	3 6 3	3 5 3	5 6 6	4 5 3	9 5 4	4 4 4	5 8 6		5 5 6	4 5 4	4 3 4		4 3 4	
5/15/10	4 3 4	5 4 3	3 4 3	0 2 2	3 4 3	9 5 4	4 4 4	5 8 6		5 5 6	4 5 4	4 3 4		4 3 4	
4/12/13	4 3 4	5 4 5	3 4 3	0 2 2	4 4 3	3 3 3	2 2 2	4 3 3		3 3 3	3 3 3	3 4 3		5 3 5	
Facial	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
Lingual															
4/12/13	6 5 4	4 3 3	4 3 4	3 2 3	4 3 4	3 3 2	3 3 3	0 2 2		3 3 3	3 4 4	3 4 3		5 4 4	
5/15/10	6 5 4	4 3 4	4 3 4	3 2 3	4 3 4	6 6 6	4 3 2	0 2 2		3 3 3	3 4 3	3 4 3		5 4 4	
5/14/10	6 5 4	4 3 4	4 3 4	3 2 3	4 3 4	6 6 6	4 3 2	0 2 2		3 3 3	3 4 3	3 4 3		5 4 4	
5/13/10	6 5 4	4 3 4	4 3 4	3 2 3	4 3 4	6 6 6	4 3 2	0 2 2		3 3 3	3 4 3	3 4 3		5 4 4	
4/30/10	6 5 4	4 3 4	4 3 4	3 2 3	4 3 4	6 6 6	4 3 2	0 2 2		3 3 3	3 4 3	3 4 3		5 4 4	
4/30/10															
5/13/10	5 6 5	5 5 5	5 5 5	4 3 4		4 3 4	4 4 4	4 4 4	4 3 4	4 3 4		4 3 4	4 3 4	5 4 5	
5/14/10	5 6 5	5 5 5	5 5 5	4 3 4		4 3 4	4 4 4	4 4 4	4 3 4	4 3 4		4 3 4	4 3 4	5 4 5	
5/15/10	5 6 5	5 5 5	5 5 5	4 3 4		4 3 4	4 4 4	4 4 4	4 3 4	4 3 4		4 3 4	4 3 4	5 4 5	
4/12/13	5 6 5	5 5 5	5 5 5	4 3 4		4 3 4	4 4 4	4 4 4	4 3 4	4 3 4		4 3 4	4 3 4	5 4 5	
Lingual															
32	31	30	29	28	27	26	25	24	23	22	21	20	19	18	17
Facial															
4/12/13	3 3 3	3 3 3	3 3 3	3 3 3		3 3 3	3 3 2	4 4 4	3 3 3	2 3 5		4 3 3	3 3 3	3 3 3	
5/15/10	5 4 5	5 4 5	4 4 4	4 3 4		4 3 4	4 3 4	4 4 4	4 4 4	2 3 7		4 3 3	3 3 3	3 3 3	
5/14/10	5 4 5	5 4 5	4 4 4	4 3 4		4 3 4	4 3 4	4 4 4	4 4 4	2 3 7		4 3 3	3 3 3	3 3 3	
5/13/10	5 4 5	5 4 5	4 4 4	4 3 4		4 3 4	4 3 4	4 4 4	4 4 4	2 3 7		4 3 3	3 3 3	3 3 3	
4/30/10	5 4 5	5 4 5	4 4 4	4 3 4		4 3 4	4 3 4	4 4 4	4 4 4	2 3 7		4 3 3	3 3 3	3 3 3	

FIGURE 5: COMPARE PREVIOUS EXAM READINGS

Print Exam Report

Print the exam report to send to the general dentist or to the periodontist. The report shows the comments, mobility, recession, facial and lingual pockets of each tooth, and a picture of each tooth with its corresponding pocket depths. If the tooth is a crown or implant, impacted or missing, that is indicated as well.

- **SELECT** *Print exam report*
- Answer "Y" to "Would you like to print a periodontal exam?"

Display Image

This is one of the most powerful features of Datacon's Perio Program. This selection displays a graphic representation of the pocket-charting screen. You can click and draw on the screen with your mouse. When you exit the display, the drawing goes away.




FIGURE 6: DISPLAY IMAGE OF PERIO CHARTING

Exam Notes

This is a free-form note area that allows you to type any notes relating to the periodontal exam. You can also customize your system to include these notes when you print the Periodontal Exam. See Miscellaneous Perio Stuff for more details.

- **Click Exam notes** and type away

Additional Exam Data

This selection allows you to edit a word processing document assigned specifically to the perio chart. It allows additional exam data specific to your office to be noted, tracked and formatted into a letter or report. For example:

- ❖ A periodontist can store specific exam findings to create letters to patients and referral sources.
- ❖ A general dentist can store a checklist to produce educational documents personalized to a patient's specific dental needs.

System Settings for Perio Charting

There are several setup options for Perio Charting that will make your perio charting life easier. You can define default comments entered for each tooth, your exam data outline, perio pocket depth, and perio report printing options.

Charting Comments

Each Perio Exam has a selection named *Comments*. These comments are user-defined and once they are set up, they are system-wide. There are 15 comments to define. Each tooth can have one comment attached.

Hint: Use the Search button to display the name of the item in *System Stuff* you are looking for.

- From the **Main Menu**, click *System Stuff>Tooth charting options, Exam data outline>Perio charting comments*
- Type the code and press **RETURN**
- Type the comment and press **RETURN**

Code	Comment
FT	Fractured Tooth
ST	Sensitive Tooth
VG	Very Guarded
I	Infected
CA	Carries
DR	Defective Restoration
BR	Bridge
EX	Recommend Extraction
FR	Frenum Attachment
NP	Not Able To Probe
CL	Clefting

Exit

FIGURE 7: PERIO CHARTING COMMENTS SETUP SCREEN

Perio Exam Report Options

Number of Copies to Print:

Suppress Display of Pocket Depths Below ☐ mm.

Moderate Pocket Depths Begin at mm.

Advanced Pocket Depths Begin at mm.

Moderate Recession Begins at mm.

Advanced Recession Begins at mm.

Text for Zero Mobility:

Exit Continue

FIGURE 8: PERIO REPORT PRINT OPTIONS

Miscellaneous Perio Stuff

Setting up Insurance Claims to Include Perio Charting

You can include the pocket charting on your electronic claims or paper claims using the QuickNotes feature.

- **Click** on *System stuff>QuickNote setup*
- **Click** the **Add** button and type “PERIO” as the QuickNote name. A note pad will then appear and you can enter the merge fields as shown in the following screen.

Name	Note
Perio	<ul style="list-style-type: none"> <perio'ul> <perio'll> <perio'lr> <perio'ur>

Exit Search Add Remove Print

FIGURE 9: PERIO QUICKNOTE SETUP

- <perio'ul> this will pull the pocket charting in the upper left quadrant
- <perio'll> this will pull the pocket charting in the lower left quadrant
- <perio'lr> this will pull the pocket charting in the lower right quadrant
- <perio'ur> this will pull the pocket charting in the upper right quadrant
- **EXIT** back to the *Main Menu*
- When posting charges, **Click Notes on form**

The screenshot shows a software window titled "Print Form for Troy Peri". Inside, there are three main sections:

- Insurance Coverage (as of 4/12/2013)**: Contains a highlighted line "Bill Prime Aetna" and a line below it that says "--- No Second Coverage ---".
- Use Treatment Plan for Pre-Determination**: Contains a highlighted line "Bill for Completed Work".
- Notes for form**: A yellow notepad area with several horizontal lines. A "QuickNote" dialog box is open over this section, with the text "perio" entered in the "QuickNote name.." field. The dialog has "Cancel", "Enter", and "Next" buttons.

At the bottom of the window are "Exit" and "Continue" buttons.

FIGURE 10: PRINTING INSURANCE FORM

- **SEARCH** and type "*PERIO*" and press **ENTER**. The note will populate with your merge fields from above.
- **Click the Remove** button to delete the quadrants you do not want to include.
- **ENTER** and finish your insurance form.