

Year-end Backup for Mac Servers

It is a good practice to make a year-end backup of your Datacon Dental System data and permanently archive it in a safe place off-site.

Which disks you need to buy:

- Buy DVDRW disks. Be sure NOT to get CDRW, DVDR, or CDR.

Put the DVDRW into your Mac:

- If you have a mouse/monitor/keyboard attached to your server, you can just eject the drive and put the disk in. Otherwise, call Datacon and we can eject the drive for you.

Where does the disk go?

See the photos below illustrating the location of the drives depending on which Mac server you have.



Mac Mini Server



Mac Pro Server

Run the backup:

- Make sure you do not have a USB drive of any kind plugged in. This includes an external backup hard drive.
- Run your backup as usual in your Datacon Dental Practice Management software.
- The DVD will eject automatically when it is completed.
- Label the DVD and store it in a safe place off-site where it will be protected from damage and heat.

Questions? Contact the staff at Datacon for help.

www.datacondental.com

Phone: (800) 773-7878

Email: info@datacondental.com